

Anti-RFC2 [R05-6U-5] Monoclonal Antibody

AGMB03896

Description

This Anti-RFC2 [R05-6U-5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03896
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	RFC2, Replication factor C subunit 2, Activator 1 40 kDa subunit (A1 40 kDa subunit), Activator 1 subunit 2, Replication factor C 40 kDa subunit (RF-C 40 kDa subunit, RFC40)
Applications:	WB ICC/IF FC
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R05-6U-5
Clonality:	Monoclonal Antibody
SwissProt ID:	P35250
Immunogen:	Peptide
Gene ID:	5982

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	RFC2
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW:39 kDa,Observed MW: 39 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	ICC/IF	1:100-1:200
	FC	1:10-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.