

Anti-Desmoglein 2/DSG2 [R01-4P-7] Monoclonal Antibody

AGMB03900

Description

This Anti-Desmoglein 2/DSG2 [R01-4P-7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03900
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	desmoglein-2, DSG2, HDGC, ARVD10, CDHF5, CMD1BB, ARVC 10, ARVC10, ARVD 10,
Applications:	WB IHC-P
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R01-4P-7
Clonality:	Monoclonal Antibody
SwissProt ID:	Q14126
Immunogen:	Peptide
Gene ID:	1829

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	DSG2
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW:122 kDa,Observed MW: 160 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:2000-1:5000
	IHC-P	1:200-1:500
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	