

# Anti-CD239/BCAM [R02-4M-1] Monoclonal Antibody

## AGMB03924

### Description

---

This Anti-CD239/BCAM [R02-4M-1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                                                                                                                                 |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB03924                                                                                                                                                                       |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)                                                                                                                             |
| <b>Synonyms:</b>      | BCAM, LU, MSK19, Basal cell adhesion molecule, Auberger B antigen, B-CAM cell surface glycoprotein, F8/G253 antigen, Lutheran antigen, Lutheran blood group glycoprotein, CD239 |
| <b>Applications:</b>  | <b>WB</b> <b>IHC-F</b> <b>IHC-P</b> <b>ICC/IF</b>                                                                                                                               |
| <b>Research Area:</b> | Cardiovascular                                                                                                                                                                  |
| <b>Form:</b>          | Liquid                                                                                                                                                                          |

### Antibody Data

---

|                      |                     |
|----------------------|---------------------|
| <b>Reactivity:</b>   | Human, Mouse, Rat   |
| <b>Clone:</b>        | R02-4M-1            |
| <b>Clonality:</b>    | Monoclonal Antibody |
| <b>SwissProt ID:</b> | P50895              |
| <b>Immunogen:</b>    | Peptide             |
| <b>Gene ID:</b>      | 4059                |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** BCAM

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW:67 kDa, Observed MW: 85kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

| <b>Antibody Dilution Ratio:</b> | <b>Application</b> | <b>Antibody Dilution Ratio</b> |
|---------------------------------|--------------------|--------------------------------|
|                                 | WB                 | 1:1000-1:5000                  |
|                                 | IHC-P              | 1:200-1:500                    |
|                                 | IF                 | 1:50-1:100                     |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.