

Anti-Cytokeratin 6A [R05-2M-5] Monoclonal Antibody

AGMB03928

Description

This Anti-Cytokeratin 6A [R05-2M-5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03928
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	KRT6A, K6A, KRT6D, Keratin, type II cytoskeletal 6A, Cytokeratin-6A (CK-6A), Cytokeratin-6D (CK-6D), Keratin-6A (K6A), Type-II keratin Kb6, Hom s 5
Applications:	WB IHC-F IHC-P ICC/IF
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R05-2M-5
Clonality:	Monoclonal Antibody
SwissProt ID:	P02538
Immunogen:	Peptide
Gene ID:	3853

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: KRT6A/KRT6B/KRT6C

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW:60 kDa, Observed MW: 60 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC-P	1:20-1:50
	ICC/IF	1:20-1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.