

# Anti-Annexin 4/ANXA4 [R01-4T-7] Monoclonal Antibody

## AGMB03931

### Description

---

This Anti-Annexin 4/ANXA4 [R01-4T-7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB03931
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	AIV, AnnexinA4, AnnexinIV, ANX 4, ANX A4, ANX4, ANXA4, Chromobindin4, EndonexinI, LipocortinIV, P32.5, P33/41, PAPII, PIG28, PP4X, Protein II, Xanx-4, ZAP36
<b>Applications:</b>	<b>WB</b>
<b>Research Area:</b>	Cardiovascular
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human, Mouse, Rat
<b>Clone:</b>	R01-4T-7
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P09525
<b>Immunogen:</b>	Peptide
<b>Gene ID:</b>	307

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** ANXA4  
**Host Species:** Rabbit  
**Isotype:** IgG  
**Purification:** Affinity Purified  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** Calculated MW:36 kDa, Observed MW: 36 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

	Application	Antibody Dilution Ratio
<b>Antibody Dilution Ratio:</b>	WB	1:1000-1:5000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.