

Anti-GITR [R09-6K-1] Monoclonal Antibody

AGMB03944

Description

This Anti-GITR [R09-6K-1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB03944
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	AITR, GITR, TNFRSF18, CD357
Applications:	FC ELISA
Research Area:	Immunology
Form:	Lyophilized

Antibody Data

Reactivity:	Human
Clone:	R09-6K-1
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9Y5U5
Immunogen:	Recombinant human GITR protein
Gene ID:	8784
Gene Name:	TNFRSF18

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at -20°C to -80°C for 12 months in lyophilized form.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	FC	1:50-1:100
	ELISA	1:10000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.