

# Anti-IL-13 Receptor alpha 1 [R05-1A-8] Monoclonal Antibody

## AGMB03964

### Description

---

This Anti-IL-13 Receptor alpha 1 [R05-1A-8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                                                                                                                                                       |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB03964                                                                                                                                                                                             |
| <b>Contents:</b>      | 100µg<br>Bradford Reagent: 1 vial (2ml)                                                                                                                                                               |
| <b>Synonyms:</b>      | IL13RA1, IL13R, IL13RA, Interleukin-13 receptor subunit alpha-1, IL-13 receptor subunit alpha-1, IL-13R subunit alpha-1, IL-13R-alpha-1, IL-13RA1, Cancer/testis antigen 19, CT19, CD antigen CD213a1 |
| <b>Applications:</b>  | <b>FC</b> <b>ELISA</b>                                                                                                                                                                                |
| <b>Research Area:</b> | Immunology                                                                                                                                                                                            |
| <b>Form:</b>          | Lyophilized                                                                                                                                                                                           |

### Antibody Data

---

|                      |                                   |
|----------------------|-----------------------------------|
| <b>Reactivity:</b>   | Human                             |
| <b>Clone:</b>        | R05-1A-8                          |
| <b>Clonality:</b>    | Monoclonal Antibody               |
| <b>SwissProt ID:</b> | P78552                            |
| <b>Immunogen:</b>    | Recombinant human IL13RA1 protein |
| <b>Gene ID:</b>      | 3597                              |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** IL13RA1

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose is added as protectants before lyophilization.

| Antibody Dilution Ratio: | Application | Antibody Dilution Ratio |
|--------------------------|-------------|-------------------------|
|                          | FC          | 1:50-1:100              |
|                          | ELISA       | 1:10000                 |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.