

Anti-alpha 1 Spectrin [R09-7X1] Monoclonal Antibody

AGMB04056

Description

This Anti-alpha 1 Spectrin [R09-7X1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04056
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	EL2, Elliptocytosis2, erythrocyte, Erythroid alpha spectrin, HS3, SPH3, SPTA, SPTA1
Applications:	WB IHC-P
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R09-7X1
Clonality:	Monoclonal Antibody
SwissProt ID:	P02549
Immunogen:	A synthesized peptide derived from human SPTA1
Gene ID:	6708

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	SPTA1
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Chromatography
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 280 kDa, Observed MW: 280 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:200-1:500
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	