

# Anti-UQCRC2 [R09-7K5] Monoclonal Antibody

## AGMB04078

### Description

---

This Anti-UQCRC2 [R09-7K5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                                                          |
|-----------------------|----------------------------------------------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB04078                                                                                                |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)                                                      |
| <b>Synonyms:</b>      | Core protein II, QCR2, mitochondrial, Ubiquinol cytochrome c reductase<br>core protein II, UQCR2, Uqcrc2 |
| <b>Applications:</b>  | <b>WB</b> <b>IHC-P</b> <b>ICC/IF</b> <b>FC</b> <b>IP</b>                                                 |
| <b>Research Area:</b> | Signal Transduction                                                                                      |
| <b>Form:</b>          | Liquid                                                                                                   |

### Antibody Data

---

|                      |                                                 |
|----------------------|-------------------------------------------------|
| <b>Reactivity:</b>   | Human, Mouse, Rat                               |
| <b>Clone:</b>        | R09-7K5                                         |
| <b>Clonality:</b>    | Monoclonal Antibody                             |
| <b>SwissProt ID:</b> | P22695                                          |
| <b>Immunogen:</b>    | A synthesized peptide derived from human UQCRC2 |
| <b>Gene ID:</b>      | 7385                                            |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

|                          |                                            |
|--------------------------|--------------------------------------------|
| <b>Gene Name:</b>        | UQCRC2                                     |
| <b>Host Species:</b>     | Rabbit                                     |
| <b>Isotype:</b>          | IgG                                        |
| <b>Purification:</b>     | Affinity Chromatography                    |
| <b>Conjugated:</b>       | Unconjugated                               |
| <b>Modification:</b>     | Unmodified                                 |
| <b>Molecular Weight:</b> | Calculated MW: 48 kDa, Observed MW: 48 kDa |

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

**Antibody Dilution Ratio:**

| Application | Antibody Dilution Ratio |
|-------------|-------------------------|
| WB          | 1:500-1:1000            |
| IHC-P       | 1:100-1:200             |
| ICC/IF      | 1:100-1:200             |
| FC          | 1:50-1:100              |
| IP          | 1:20-1:50               |

**Protein Quantification (Optional):**

To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.