

Anti-LAMP2 [R01-2P3] Monoclonal Antibody

AGMB04084

Description

This Anti-LAMP2 [R01-2P3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04084
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	LAMP2, Lysosome-associated membrane glycoprotein 2, LAMP-2, Lysosome-associated membrane protein 2, CD107 antigen-like family member B, CD107b, CD107 antigen-like family member B, CD107b, LAMP 2,
Applications:	WB IHC-P IP
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R01-2P3
Clonality:	Monoclonal Antibody
SwissProt ID:	P13473
Immunogen:	A synthesized peptide derived from human LAMP2a

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 3920

Gene Name: LAMP2

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 45 kDa, Observed MW: 120 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:100-1:200
	IP	1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.