

Anti-IKK alpha/beta [R01-3A4] Monoclonal Antibody

AGMB04095

Description

This Anti-IKK alpha/beta [R01-3A4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04095
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CHUK, IKKA, TCF16, Inhibitor of nuclear factor kappa-B kinase subunit alpha, I-kappa-B kinase alpha, IKK-A, IKK-alpha, IκBKA, IκappaB kinase, Conserved helix-loop-helix ubiquitous kinase, I-kappa-B kinase 1, IKK1, Nuclear factor NF-kappa-B
Applications:	WB ICC/IF IP
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R01-3A4
Clonality:	Monoclonal Antibody
SwissProt ID:	O15111/O14920
Immunogen:	A synthesized peptide derived from human IKK alpha+beta

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 1147/3551

Gene Name: CHUK/IKBKB

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 85 kDa, Observed MW: 75,87 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:

Application	Antibody Dilution Ratio
WB	1:1000-1:2000
IF	1:50-1:200
IP	1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.