

# Anti-beta 2 Adrenergic Receptor [R06-7A7] Monoclonal Antibody

## AGMB04106

### Description

---

This Anti-beta 2 Adrenergic Receptor [R06-7A7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                                             |
|-----------------------|---------------------------------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB04106                                                                                   |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)                                         |
| <b>Synonyms:</b>      | ADRB2, ADRB2R, B2AR, Beta-2 adrenergic receptor, Beta-2 adrenoreceptor, Beta-2 adrenoceptor |
| <b>Applications:</b>  | <b>WB</b> <b>IHC-P</b> <b>IP</b>                                                            |
| <b>Research Area:</b> | Signal Transduction                                                                         |
| <b>Form:</b>          | Liquid                                                                                      |

### Antibody Data

---

|                      |                                                                     |
|----------------------|---------------------------------------------------------------------|
| <b>Reactivity:</b>   | Human, Mouse, Rat                                                   |
| <b>Clone:</b>        | R06-7A7                                                             |
| <b>Clonality:</b>    | Monoclonal Antibody                                                 |
| <b>SwissProt ID:</b> | P07550                                                              |
| <b>Immunogen:</b>    | A synthesized peptide derived from human beta 2 Adrenergic Receptor |
| <b>Gene ID:</b>      | 154                                                                 |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** ADRB2

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 46 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

| <b>Antibody Dilution Ratio:</b> | <b>Application</b> | <b>Antibody Dilution Ratio</b> |
|---------------------------------|--------------------|--------------------------------|
|                                 | WB                 | 1:1000-1:2000                  |
|                                 | IHC                | 1:100-1:200                    |
|                                 | IP                 | 1:30                           |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.