

Anti-CXCL11 [R07-5Q1] Monoclonal Antibody

AGMB04112

Description

This Anti-CXCL11 [R07-5Q1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04112
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Beta-R1, CXC11, CXCL11, H174, I TAC, IP9, ITAC, SCYB11, SCYB9B, Interferon gamma-inducible protein 9 (IP-9), Interferon-inducible T-cell alpha chemoattractant (I-TAC), Small-inducible cytokine B11
Applications:	WB ICC/IF
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R07-5Q1
Clonality:	Monoclonal Antibody
SwissProt ID:	O14625
Immunogen:	A synthesized peptide derived from human CXCL11
Gene ID:	6373

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CXCL11

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 10 kDa, Observed MW: 15 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.