

# Anti-KLF4 [R03-5H7] Monoclonal Antibody

AGMB04150

## Description

---

This Anti-KLF4 [R03-5H7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB04150
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	Endothelial Kruppel like zinc finger protein, Epithelial zinc finger protein EZF, EZF, GKLF, Gut-enriched krueppel-like factor, KLF, KLF4, KLF4_HUMAN, Krueppel-like factor 4, Kruppel like factor 4, Epithelial zinc finger protein EZF, Gut enriched Krueppel like factor, Kruppel like factor 4 (gut).
<b>Applications:</b>	<b>WB</b> <b>IHC-P</b> <b>ICC/IF</b>
<b>Research Area:</b>	Epigenetics
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Human, Mouse, Rat
<b>Clone:</b>	R03-5H7
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	O43474
<b>Immunogen:</b>	A synthesized peptide derived from human KLF4

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene ID:** 9314  
**Gene Name:** KLF4  
**Host Species:** Rabbit  
**Isotype:** IgG  
**Purification:** Affinity Chromatography  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** Calculated MW: 55 kDa, Observed MW: 55 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

**Antibody Dilution Ratio:**

Application	Antibody Dilution Ratio
WB	1:1000-1:2000
IHC	1:100-1:200
IF	1:50-1:200

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.