

# Anti-P Glycoprotein [R05-9K3] Monoclonal Antibody

## AGMB04154

### Description

---

This Anti-P Glycoprotein [R05-9K3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB04154
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	p-pg, ABCB1, MDR1, PGY1, Multidrug resistance protein 1, ATP-binding cassette sub-family B member 1, P-glycoprotein 1, CD antigen CD243, ABC20, ABCB1, ATP binding cassette, sub family B (MDR/TAP), member 1,
<b>Applications:</b>	<b>WB</b> <b>IHC-P</b>
<b>Research Area:</b>	Signal Transduction
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human, Mouse, Rat
<b>Clone:</b>	R05-9K3
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P08183
<b>Immunogen:</b>	A synthesized peptide derived from human P Glycoprotein
<b>Gene ID:</b>	5243

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** ABCB1  
**Host Species:** Rabbit  
**Isotype:** IgG  
**Purification:** Affinity Chromatography  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** Calculated MW: 141 kDa, Observed MW: 180 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

	Application	Antibody Dilution Ratio
<b>Antibody Dilution Ratio:</b>	WB	1:1000-1:2000
	IHC	1:50-1:100

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.