

Anti-RIP3 (6B4) [6B4-8G2-10A2] Monoclonal Antibody

AGMB04170

Description

This Anti-RIP3 (6B4) [6B4-8G2-10A2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04170
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	RIP3, RIP like protein kinase 3, RIP 3, Receptor-interacting serine/threonine-protein kinase 3, Receptor-interacting protein 3
Applications:	WB ELISA
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	6B4-8G2-10A2
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9Y572
Immunogen:	This RIPK3 antibody is generated from a mouse immunized with a recombinant protein from the human region of human RIPK3.
Gene ID:	11035

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: RIPK3

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 57 kDa, Observed MW: 57 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 0.09% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	ELISA	1:10000
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	