

Anti-Lin28A (2C1) [2C1-F9-A2] Monoclonal Antibody

AGMB04180

Description

This Anti-Lin28A (2C1) [2C1-F9-A2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04180
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Protein lin-28 homolog A, LIN28A, CSDD1, LIN28, ZCCHC1, CSDD1, FLJ12457, LIN 28,
Applications:	WB
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse
Clone:	2C1-F9-A2
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9H9Z2
Immunogen:	Purified recombinant human LIN28A protein fragments expressed in E.coli
Gene ID:	79727

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: LIN28A
Host Species: Mouse
Isotype: IgG2a
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 23 kDa, Observed MW: 26 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.