

# Anti-LC3A/B (3E9) [3E9-E5-C9] Monoclonal Antibody

## AGMB04183

### Description

---

This Anti-LC3A/B (3E9) [3E9-E5-C9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB04183
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	LC3, LC3A, ATG8E, MAP1ALC3, MAP1BLC3, MAP1LC3A, LC3B, ATG8F, MAP1LC3B-a, MAP1A/1BLC3, MAP1LC3B
<b>Applications:</b>	<b>WB</b>
<b>Research Area:</b>	Signal Transduction
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human, Rat
<b>Clone:</b>	3E9-E5-C9
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	Q9H492/Q9GZQ8
<b>Immunogen:</b>	Synthetic peptide corresponding to human LC3B protein
<b>Gene ID:</b>	84557/81631

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** MAP1LC3A/MAP1LC3B

**Host Species:** Mouse

**Isotype:** IgG2b

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 14 kDa, Observed MW: 14,16 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

<b>Antibody Dilution Ratio:</b>	<b>Application</b>	<b>Antibody Dilution Ratio</b>
	WB	1:500-1:1000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.