

Anti-LPP (6F6) [6F6-C1-B6] Monoclonal Antibody

AGMB04199

Description

This Anti-LPP (6F6) [6F6-C1-B6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04199
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	LPP, Lipoma-preferred partner, LIM domain-containing preferred translocation partner in lipoma
Applications:	WB IHC-F IHC-P ICC/IF IP
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat, Hamster, Monkey
Clone:	6F6-C1-B6
Clonality:	Monoclonal Antibody
SwissProt ID:	Q93052
Immunogen:	Purified recombinant fragment of human LPP expressed in E. Coli.
Gene ID:	4026

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	LPP
Host Species:	Mouse
Isotype:	IgG1
Purification:	Ascitic Fluid
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 66 kDa, Observed MW: 66 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Ascitic fluid containing 0.03% sodium azide.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:50-1:100
	IF	1:50-1:200
	IP	1:20
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	