

Anti-SMYD1 (2E11) [2E11-B6-C7] Monoclonal Antibody

AGMB04205

Description

This Anti-SMYD1 (2E11) [2E11-B6-C7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04205
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	BOP, KMT3D, ZMYND18, ZMYND22, CD8 beta opposite, CD8b opposite,
Applications:	
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Rat
Clone:	2E11-B6-C7
Clonality:	Monoclonal Antibody
SwissProt ID:	Q8NB12
Immunogen:	Purified recombinant human SMYD1 protein fragments expressed in E.coli.
Gene ID:	150572

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: SMYD1
Host Species: Mouse
Isotype: IgG2b
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 57 kDa, Observed MW: 57 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.