

Anti-SAFB1 (2E8) [2E8-E2-G6] Monoclonal Antibody

AGMB04227

Description

This Anti-SAFB1 (2E8) [2E8-E2-G6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04227
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	SAFB, HAP, HET, SAFB1, Scaffold attachment factor B1, SAF-B, SAF-B1, HSP27 estrogen response element-TATA box-binding protein, HSP27 ERE-TATA-binding protein
Applications:	WB ICC/IF CHIP
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	2E8-E2-G6
Clonality:	Monoclonal Antibody
SwissProt ID:	Q15424
Immunogen:	Purified recombinant human SAFB1 protein fragments expressed in E.coli.
Gene ID:	6294

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: SAFB
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 103 kDa, Observed MW: 130 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IF	1:50-1:200
	ChIP	1:20

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.