

Anti-cIAP2 (1F12) [1F12-E5-G12] Monoclonal Antibody

AGMB04238

Description

This Anti-cIAP2 (1F12) [1F12-E5-G12] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04238
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	cIAP2, AIP1, API2, BIRC3, CIAP2, HAIP1, IAP homolog C, RNF49, MIHC, MALT2, Apoptosis inhibitor 2, IAP2
Applications:	WB
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Monkey
Clone:	1F12-E5-G12
Clonality:	Monoclonal Antibody
SwissProt ID:	Q13489
Immunogen:	Purified recombinant human c-IAP2 protein fragments expressed in E.coli
Gene ID:	330

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: BIRC3
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 68 kDa, Observed MW: 72 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.