

Anti-Cytokeratin 17 (8C6) [8C6-4H2-5C1] Monoclonal Antibody

AGMB04253

Description

This Anti-Cytokeratin 17 (8C6) [8C6-4H2-5C1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04253
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Keratin, type I cytoskeletal 17, Cytokeratin-17, CK-17, Keratin-17, K17
Applications:	WB IHC-P ICC/IF IP
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	8C6-4H2-5C1
Clonality:	Monoclonal Antibody
SwissProt ID:	Q04695
Immunogen:	Synthetic Peptide of CK17
Gene ID:	3872
Gene Name:	KRT17

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Mouse
Isotype:	IgG1
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 48 kDa, Observed MW: 48 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:50-1:100
	IF	1:50-1:200
	IP	1:20

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.