

# Anti-NF-KB p65 (3D2) [3D2-4E9-7A8] Monoclonal Antibody

## AGMB04255

### Description

---

This Anti-NF-KB p65 (3D2) [3D2-4E9-7A8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB04255
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	NFKB3, RELA, TF65, Transcription factor p65, p65, NFkB
<b>Applications:</b>	<b>WB</b> <b>IHC-F</b> <b>IHC-P</b> <b>ICC/IF</b> <b>IP</b>
<b>Research Area:</b>	Cell Biology
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human, Mouse, Rat
<b>Clone:</b>	3D2-4E9-7A8
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	Q04206
<b>Immunogen:</b>	Recombinant Protein of Transcription factor p65
<b>Gene ID:</b>	5970
<b>Gene Name:</b>	RELA

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** Mouse

**Isotype:** IgG1

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 60 kDa, Observed MW: 65 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 1xPBS(pH 7.4), 150mM NaCl, 50% Glycerol, 0.02% Sodium azide and 0.05% BSA

<b>Antibody Dilution Ratio:</b>	<b>Application</b>	<b>Antibody Dilution Ratio</b>
	WB	1:500-1:1000
	IHC	1:50-1:100
	IF	1:50-1:200
	IP	1:20

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.