

Anti-GRB2 (5G3) [5G3-E7-F11] Monoclonal Antibody

AGMB04272

Description

This Anti-GRB2 (5G3) [5G3-E7-F11] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04272
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	GRB2, ASH, Growth factor receptor-bound protein 2, Adapter protein GRB2, Protein Ash, SH2/SH3 adapter GRB2
Applications:	
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	5G3-E7-F11
Clonality:	Monoclonal Antibody
SwissProt ID:	P62993
Immunogen:	Recombinant protein of human GRB2.
Gene ID:	2885

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: GRB2
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 25 kDa, Observed MW: 25 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.