

Anti-RPS6 (7B10) [7B10-E11-E11] Monoclonal Antibody

AGMB04273

Description

This Anti-RPS6 (7B10) [7B10-E11-E11] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04273
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	RPS6, OK/SW-cl.2, 40S ribosomal protein S6, Phosphoprotein NP33
Applications:	WB ICC/IF IP
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat, Monkey
Clone:	7B10-E11-E11
Clonality:	Monoclonal Antibody
SwissProt ID:	P62753
Immunogen:	Recombinant full length protein corresponding to human S6 Ribosomal Protein.
Gene ID:	6194

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: RPS6
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 29 kDa, Observed MW: 32 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IF	1:50-1:200
	IP	1:20

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.