

Anti-Aldehyde dehydrogenase 10 [3F7-H9-E12] Monoclonal Antibody

AGMB04293

Description

This Anti-Aldehyde dehydrogenase 10 [3F7-H9-E12] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04293
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ALDH3A2, ALDH10, FALDH, Fatty aldehyde dehydrogenase, Aldehyde dehydrogenase 10, Aldehyde dehydrogenase family 3 member A2, Microsomal aldehyde dehydrogenase
Applications:	WB
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	3F7-H9-E12
Clonality:	Monoclonal Antibody
SwissProt ID:	P51648
Immunogen:	Purified recombinant human Aldehyde dehydrogenase 10 protein fragments expressed in E.coli.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 224
Gene Name: ALDH3A2
Host Species: Mouse
Isotype: IgG2a
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 55 kDa, Observed MW: 55 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.