

Anti-PI3 Kinase p85 alpha (1C8) [1C8-2H4-9D9] Monoclonal Antibody

AGMB04347

Description

This Anti-PI3 Kinase p85 alpha (1C8) [1C8-2H4-9D9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04347
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	PIK3R1, GRB1, p85, AGM7, GRB1, IMD36, p85alpha, p85-ALPHA, Phosphatidylinositol 3-kinase regulatory subunit alpha, PI3-kinase regulatory subunit alpha, PI3K regulatory subunit alpha, PtdIns-3-kinase regulatory subunit alpha,
Applications:	WB IHC-P
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	1C8-2H4-9D9
Clonality:	Monoclonal Antibody
SwissProt ID:	P27986
Immunogen:	Purified recombinant protein expressed in E.coli.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	5295
Gene Name:	PIK3R1
Host Species:	Mouse
Isotype:	IgG1
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 84 kDa, Observed MW: 84 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:50-1:100
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	