

# Anti-S100 beta (5G3) [5G3-C11] Monoclonal Antibody

## AGMB04365

### Description

---

This Anti-S100 beta (5G3) [5G3-C11] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB04365
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	NEF, S100, S100-B, S100beta, S100B
<b>Applications:</b>	<b>WB</b>
<b>Research Area:</b>	Cell Biology
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Rat
<b>Clone:</b>	5G3-C11
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P04271
<b>Immunogen:</b>	Human S100B synthetic peptide conjugated to KLH.
<b>Gene ID:</b>	6285
<b>Gene Name:</b>	S100B

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** Mouse

**Isotype:** IgG2b

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 11 kDa, Observed MW: 11 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.