

Anti-Human Serum Albumin (1A8) [1A8-C5-B4] Monoclonal Antibody

AGMB04372

Description

This Anti-Human Serum Albumin (1A8) [1A8-C5-B4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04372
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ALB, ALBU_HUMAN, Albumin (32 AA), Albumin (AA 34), Albumin, Analbuminemia, Bisalbuminemia, Cell growth inhibiting protein 42, DKFZp779N1935, Dysalbuminemic hyperthyroxinemia, Growth inhibiting protein 20, HSA, Hyperthyroxinemia dysalbuminemic, Serum albumin.
Applications:	WB ELISA
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	1A8-C5-B4
Clonality:	Monoclonal Antibody
SwissProt ID:	P02768
Immunogen:	Purified Human serum albumin.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	213
Gene Name:	ALB
Host Species:	Mouse
Isotype:	IgG1
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 69 kDa, Observed MW: 69 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	ELISA	1:10000
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	