

Anti-Cystatin C (6A4) [6A4-8C9-4B1] Monoclonal Antibody

AGMB04395

Description

This Anti-Cystatin C (6A4) [6A4-8C9-4B1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04395
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CST3, Cystatin-C, Cystatin-3, Gamma-trace, Neuroendocrine basic polypeptide, Post-gamma-globulin, AD 8, AD8, Amyloid angiopathy and cerebral hemorrhage,
Applications:	WB IHC-P ELISA
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	6A4-8C9-4B1
Clonality:	Monoclonal Antibody
SwissProt ID:	P01034
Immunogen:	Purified recombinant protein expressed in E.coli.
Gene ID:	1471

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	CST3
Host Species:	Mouse
Isotype:	IgG1
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 16 kDa, Observed MW: 16 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:50-1:100
	ELISA	1:10000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.