

Anti-ITGA2B/CD41 (9G4) [9G4-3F4-10B8] Monoclonal Antibody

AGMB04419

Description

This Anti-ITGA2B/CD41 (9G4) [9G4-3F4-10B8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04419
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	GT, GT1, GTA, CD41, GP2B, HPA3, CD41B, GPIIb, BDPLT2, BDPLT16, PPP1R93, CD41/ITGA2B
Applications:	IHC-P
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	9G4-3F4-10B8
Clonality:	Monoclonal Antibody
SwissProt ID:	P08514
Immunogen:	Synthetic Peptide of CD41
Gene ID:	3674

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: ITGA2B
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.