

Anti-CD2 (3F10) [3F10-7F8-3A1] Monoclonal Antibody

AGMB04422

Description

This Anti-CD2 (3F10) [3F10-7F8-3A1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04422
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CD2, SRBC, T-cell surface antigen CD2, Erythrocyte receptor, LFA-2, LFA-3 receptor, Rosette receptor, T-cell surface antigen T11/Leu-5, CD2
Applications:	IHC-P
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	3F10-7F8-3A1
Clonality:	Monoclonal Antibody
SwissProt ID:	P06729
Immunogen:	Synthetic peptide conjugated to KLH.
Gene ID:	914

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CD2
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.