

Anti-CD5 (10F9) [10F9-5E2-3D3] Monoclonal Antibody

AGMB04423

Description

This Anti-CD5 (10F9) [10F9-5E2-3D3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04423
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CD5, LEU1, T-cell surface glycoprotein CD5, Lymphocyte antigen T1/Leu-1, CD antigen CD5
Applications:	IHC-P ICC/IF
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	10F9-5E2-3D3
Clonality:	Monoclonal Antibody
SwissProt ID:	P06127
Immunogen:	Synthetic peptide conjugated to KLH.
Gene ID:	921

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CD5
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

	Application	Antibody Dilution Ratio
Antibody Dilution Ratio:	IHC	1:50-1:100
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.