

# Anti-epsilon Tubulin (4C4) [4C4-4H9-8C8] Monoclonal Antibody

## AGMB04432

### Description

---

This Anti-epsilon Tubulin (4C4) [4C4-4H9-8C8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                     |
|-----------------------|-----------------------------------------------------|
| <b>SKU:</b>           | AGMB04432                                           |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml) |
| <b>Synonyms:</b>      | FLJ22589, TUBE1, tubulin epsilon 1, ε tubulin       |
| <b>Applications:</b>  | <b>IHC-P</b>                                        |
| <b>Research Area:</b> | Signal Transduction                                 |
| <b>Form:</b>          | Liquid                                              |

### Antibody Data

---

|                      |                                      |
|----------------------|--------------------------------------|
| <b>Reactivity:</b>   | Human, Mouse, Rat                    |
| <b>Clone:</b>        | 4C4-4H9-8C8                          |
| <b>Clonality:</b>    | Monoclonal Antibody                  |
| <b>SwissProt ID:</b> | Q9UJT0                               |
| <b>Immunogen:</b>    | Synthetic peptide conjugated to KLH. |
| <b>Gene ID:</b>      | 51175                                |
| <b>Gene Name:</b>    | TUBE1                                |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** Mouse

**Isotype:** IgG1

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

| <b>Antibody Dilution Ratio:</b> | <b>Application</b> | <b>Antibody Dilution Ratio</b> |
|---------------------------------|--------------------|--------------------------------|
|                                 | IHC                | 1:50-1:100                     |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.