

## Anti-HP1 gamma (9E5) [9E5-1E3-3A8] Monoclonal Antibody

AGMB04465

### Description

---

This Anti-HP1 gamma (9E5) [9E5-1E3-3A8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                                                                 |
|-----------------------|-----------------------------------------------------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB04465                                                                                                       |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)                                                             |
| <b>Synonyms:</b>      | CBX3, Chromobox protein homolog 3, HECH, Heterochromatin protein 1 homolog gamma, HP1 gamma, Modifier 2 protein |
| <b>Applications:</b>  | <b>WB</b> <b>IHC-P</b>                                                                                          |
| <b>Research Area:</b> | Epigenetics                                                                                                     |
| <b>Form:</b>          | Liquid                                                                                                          |

### Antibody Data

---

|                      |                                                   |
|----------------------|---------------------------------------------------|
| <b>Reactivity:</b>   | Human, Mouse, Rat                                 |
| <b>Clone:</b>        | 9E5-1E3-3A8                                       |
| <b>Clonality:</b>    | Monoclonal Antibody                               |
| <b>SwissProt ID:</b> | Q13185                                            |
| <b>Immunogen:</b>    | Purified recombinant protein expressed in E.coli. |
| <b>Gene ID:</b>      | 11335                                             |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** CBX3  
**Host Species:** Mouse  
**Isotype:** IgG1  
**Purification:** Affinity Purified  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** Calculated MW: 21 kDa, Observed MW: 24 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

| Antibody Dilution Ratio: | Application | Antibody Dilution Ratio |
|--------------------------|-------------|-------------------------|
|                          | WB          | 1:500-1:1000            |
|                          | IHC         | 1:50-1:100              |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.