

# Anti-FABP3 (3E8B10) [3E8-B10] Monoclonal Antibody

AGMB04478

## Description

---

This Anti-FABP3 (3E8B10) [3E8-B10] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB04478
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	422 protein, Cardiac Fatty Acid Binding Protein, FABP 11, FABP 3, FABP11, FABP3, FABPH_HUMAN, fatty acid binding protein 11, Fatty acid binding protein 3, Fatty acid binding protein 3 muscle and heart, Fatty acid binding protein 3 muscle and heart mammary derived growth inhibitor, Fatty acid binding protein 3 muscle, Fatty acid binding protein 3
<b>Applications:</b>	<b>ELISA</b>
<b>Research Area:</b>	Cardiovascular
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	3E8-B10
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P05413
<b>Immunogen:</b>	Purified recombinant H-FABP protein expressed in E.coli

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene ID:** 2170  
**Gene Name:** FABP3  
**Host Species:** Mouse  
**Isotype:** IgG2b  
**Purification:** Affinity Purified  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	ELISA	1:10000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.