

Anti-Leptin (5C11) [5C11-G6-B12] Monoclonal Antibody

AGMB04481

Description

This Anti-Leptin (5C11) [5C11-G6-B12] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04481
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	LEP, Leptin Murine Obesity Homolog, Leptin Precursor Obesity Factor, OB, Obese Protein, Obesity, Obesity factor, Obesity homolog mouse, Obesity Murine Homolog Leptin, OBS, LEP_HUMAN.
Applications:	ELISA
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	5C11-G6-B12
Clonality:	Monoclonal Antibody
SwissProt ID:	P41159
Immunogen:	Purified recombinant full length of human Leptin protein expressed in E.coli.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 3952

Gene Name: LEP

Host Species: Mouse

Isotype: IgG2a

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	ELISA	1:10000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.