

Anti-RBP4 (7C5) [7C5-F6] Monoclonal Antibody

AGMB04487

Description

This Anti-RBP4 (7C5) [7C5-F6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04487
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	OTTHUMP00000020114, OTTHUMP00000020115, OTTHUMP00000020116, Plasma retinol binding protein 4, Plasma retinol-binding protein, Plasma retinol-binding protein(1-176), PRBP, PRO2222, RBP, RBP4, RDCCAS, RET4_HUMAN, Retinol binding protein 4, Retinol binding protein 4 interstitial, Retinol binding protein 4 plasma.
Applications:	ELISA
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	7C5-F6
Clonality:	Monoclonal Antibody
SwissProt ID:	P02753
Immunogen:	Purified recombinant human RBP4 protein fragments expressed in E.coli.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 5950
Gene Name: RBP4
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	ELISA	1:10000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.