

# Anti-TBL1X (4H2) [4H2-D5-E9] Monoclonal Antibody

## AGMB04502

### Description

---

This Anti-TBL1X (4H2) [4H2-D5-E9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB04502
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	TBL1X, TBL1, F-box-like/WD repeat-containing protein TBL1X, SMAP55, Transducin beta-like protein 1X, Transducin-beta-like protein 1, X-linked
<b>Applications:</b>	<b>WB</b> <b>ICC/IF</b>
<b>Research Area:</b>	Cell Biology
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	4H2-D5-E9
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	O60907
<b>Immunogen:</b>	Purified recombinant human TBL1 protein fragments expressed in E.coli.
<b>Gene ID:</b>	6907

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** TBL1X  
**Host Species:** Mouse  
**Isotype:** IgG2b  
**Purification:** Affinity Purified  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** Calculated MW: 62 kDa, Observed MW: 58 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

	Application	Antibody Dilution Ratio
<b>Antibody Dilution Ratio:</b>	WB	1:500-1:1000
	IF	1:50-1:200

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.