

Anti-EEF2/Elongation factor 2 (4B3) [4B3-G7-H5] Monoclonal Antibody AGMB04515

Description

This Anti-EEF2/Elongation factor 2 (4B3) [4B3-G7-H5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04515
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	EEF2, EF2, Elongation factor 2, EF-2
Applications:	WB ICC/IF CHIP
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	4B3-G7-H5
Clonality:	Monoclonal Antibody
SwissProt ID:	P13639
Immunogen:	Purified recombinant human eEF2 protein fragments expressed in E.coli.
Gene ID:	1938
Gene Name:	EEF2

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Mouse

Isotype: IgG2b

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 95 kDa, Observed MW: 95 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IF	1:50-1:200
	ChIP	1:20

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.