

# Anti-HDAC3 (3G3) [3G3-H6-H10] Monoclonal Antibody

## AGMB04528

### Description

---

This Anti-HDAC3 (3G3) [3G3-H6-H10] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                                     |
|-----------------------|-------------------------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB04528                                                                           |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)                                 |
| <b>Synonyms:</b>      | Histone deacetylase 3, HD3, Protein deacetylase HDAC3, HDAC3                        |
| <b>Applications:</b>  |  |
| <b>Research Area:</b> | Epigenetics                                                                         |
| <b>Form:</b>          | Liquid                                                                              |

### Antibody Data

---

|                      |                                                                         |
|----------------------|-------------------------------------------------------------------------|
| <b>Reactivity:</b>   | Human, Mouse, Rat                                                       |
| <b>Clone:</b>        | 3G3-H6-H10                                                              |
| <b>Clonality:</b>    | Monoclonal Antibody                                                     |
| <b>SwissProt ID:</b> | O15379                                                                  |
| <b>Immunogen:</b>    | Purified recombinant human HDAC3 protein fragments expressed in E.coli. |
| <b>Gene ID:</b>      | 8841                                                                    |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** HDAC3  
**Host Species:** Mouse  
**Isotype:** IgG1  
**Purification:** Affinity Purified  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** Calculated MW: 49 kDa, Observed MW: 49 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

| Antibody Dilution Ratio: | Application | Antibody Dilution Ratio |
|--------------------------|-------------|-------------------------|
|                          | WB          | 1:500-1:1000            |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.