

Anti-SMYD4 (7C8) [7C8-B11-F10] Monoclonal Antibody

AGMB04538

Description

This Anti-SMYD4 (7C8) [7C8-B11-F10] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04538
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ZMYND21, SMYD 4, SET and MYND domain-containing protein 4
Applications:	WB
Research Area:	Cancer
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	7C8-B11-F10
Clonality:	Monoclonal Antibody
SwissProt ID:	Q8IYR2
Immunogen:	Purified recombinant human SMYD4 protein fragments expressed in E.coli.
Gene ID:	114826

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: SMYD4
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 89 kDa, Observed MW: 89 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.