

Anti-S100A10 (6F4) [6F4-E6-D5-C10] Monoclonal Antibody

AGMB04561

Description

This Anti-S100A10 (6F4) [6F4-E6-D5-C10] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04561
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	S100A10, 42C, ANX2L, ANX2LG, CAL1L, CLP11, Ca[1], GP11, P11, p10, AA409961, AL024248,
Applications:	WB ICC/IF IP
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	6F4-E6-D5-C10
Clonality:	Monoclonal Antibody
SwissProt ID:	P08207
Immunogen:	Purified recombinant mouse S100A10 protein fragments expressed in E.coli.
Gene ID:	20194

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: S100a10

Host Species: Mouse

Isotype: IgG2a

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 11 kDa, Observed MW: 11 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IF	1:50-1:200
	IP	1:20

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.