

# Anti-Cyclophilin B [6H1-10C4-1A2] Monoclonal Antibody

## AGMB04591

### Description

---

This Anti-Cyclophilin B [6H1-10C4-1A2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB04591
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	PPIB, CYPB, Peptidyl-prolyl cis-trans isomerase B, PPlase B, CYP-S1, Cyclophilin B, Rotamase B, S-cyclophilin, SCYLP
<b>Applications:</b>	<b>WB</b>
<b>Research Area:</b>	Immunology
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human, Mouse, Rat
<b>Clone:</b>	6H1-10C4-1A2
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P23284
<b>Immunogen:</b>	Synthetic peptide conjugated to KLH.
<b>Gene ID:</b>	5479

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** PPIB  
**Host Species:** Mouse  
**Isotype:** IgG1  
**Purification:** Affinity Purified  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** Calculated MW: 24 kDa, Observed MW: 24 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

<b>Antibody Dilution Ratio:</b>	<b>Application</b>	<b>Antibody Dilution Ratio</b>
	WB	1:500-1:1000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.