

Anti-NSE [1H8-B4] Monoclonal Antibody

AGMB04610

Description

This Anti-NSE [1H8-B4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04610
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ENO2, Gamma-enolase, 2-phospho-D-glycerate hydro-lyase, Enolase 2, Neural enolase, Neuron-specific enolase, NSE
Applications:	WB ELISA
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	1H8-B4
Clonality:	Monoclonal Antibody
SwissProt ID:	P09104
Immunogen:	Synthetic peptide derived from human NSE.
Gene ID:	2026

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: ENO2
Host Species: Mouse
Isotype: IgG1
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 47 kDa, Observed MW: 47 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

	Application	Antibody Dilution Ratio
Antibody Dilution Ratio:	WB	1:500-1:1000
	ELISA	1:10000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.