

Anti-Cardiac Troponin I (2D10) [2D10-C5] Monoclonal Antibody

AGMB04612

Description

This Anti-Cardiac Troponin I (2D10) [2D10-C5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04612
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	TNNI3, TNNC1, Troponin I, cardiac muscle, Cardiac troponin I
Applications:	WB ELISA
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	2D10-C5
Clonality:	Monoclonal Antibody
SwissProt ID:	P19429
Immunogen:	Synthetic peptides corresponding to the sequence of human Cardiac Troponin I.
Gene ID:	7137

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	TNNI3
Host Species:	Mouse
Isotype:	IgG2a
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 24 kDa, Observed MW: 28 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	ELISA	1:10000
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	