

Anti-Lgi4 Rat [11H12] Monoclonal Antibody

AGMB04656

Description

This Anti-Lgi4 Rat [11H12] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04656
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	LGI1-like protein 3, Leucine-rich glioma-inactivated protein 4
Applications:	WB IHC-P
Research Area:	-
Form:	Liquid

Antibody Data

Reactivity:	Mouse
Clone:	11H12
Clonality:	Monoclonal Antibody
SwissProt ID:	Q8K1S1
Immunogen:	Mouse Lgi4 peptide coupled to carrier protein
Gene ID:	243914
Gene Name:	Lgi4

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Rat
Isotype:	IgG1
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 59 kDa, Observed MW: 59 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 0.09% sodium azide as a preservative.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	