

Anti-Phospho-EB3 (Ser162) [4B7] Monoclonal Antibody

AGMB04657

Description

This Anti-Phospho-EB3 (Ser162) [4B7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04657
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	EB3, EBF3, MAPRE3, RP/EB Family, RP3, APC binding protein, EB 3, EB1 protein family member 3,
Applications:	WB IHC-P
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	4B7
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9UPY8
Immunogen:	S162 phosphorylated EB3 peptide conjugated to carrier protein
Gene ID:	22924

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	MAPRE3
Host Species:	Mouse
Isotype:	IgG2b
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 32 kDa, Observed MW: 32 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 0.09% sodium azide as a preservative.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:50-1:100
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	