

Anti-HIV 1 p24 (1F8) [1F8-B5] Monoclonal Antibody

AGMB04668

Description

This Anti-HIV 1 p24 (1F8) [1F8-B5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04668
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Human Immunodeficiency Virus Type 1 (HIV-1) gag p24
Applications:	ELISA
Research Area:	Microbiology
Form:	Liquid

Antibody Data

Reactivity:	Virus
Clone:	1F8-B5
Clonality:	Monoclonal Antibody
SwissProt ID:	P04591
Immunogen:	Recombinant full-length HIV-1 p24 protein.
Gene ID:	155030
Gene Name:	gag

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	ELISA	1:10000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.